

NAT2 Antibody



PACO10740

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NAT2

Uniprot

P11245

Synonyms:

N-acetyltransferase 2 (arylamine N-acetyltransferase); NAT2; AAC2; PNAT ;

Immunogen:

Human NAT2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A