

# NCAPD3 Antibody



PACO10752

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NCAPD3

**Uniprot**

P42695

**Synonyms:**

non-SMC condensin II complex, subunit D3; NCAPD3; CAP-D3; FLJ42888; KIAA0056; MGC104671; hCAP-D3; hHCP-6; hcp-6 ;

**Immunogen:**

Human NCAPD3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A