

PACO10754

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NCF4

Uniprot

Q15080

Synonyms:

neutrophil cytosolic factor 4, 40kDa; NCF4; MGC3810; NCF; P40PHOX; SH3PXD4 ;

Immunogen:

Human NCF4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A