

NCKAP1 Antibody



PACO10757

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NCKAP1

Uniprot

Q9Y2A7

Synonyms:

NCK-associated protein 1; NCKAP1; FLJ11291; HEM2; KIAA0587; MGC8981; NAP1; NAP125 ;

Immunogen:

Human NCKAP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A