

PACO10765

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NCOR1

Uniprot

O75376

Synonyms:

nuclear receptor co-repressor 1; NCOR1; KIAA1047; MGC104216; N-CoR; TRAC1; hCIT529110; hN-CoR ;

Immunogen:

Human NCOR1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A