

NDOR1 Antibody



PACO10770

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NDOR1

Uniprot

Q9UHB4

Synonyms:

NADPH dependent diflavin oxidoreductase 1; NDOR1; MGC138148; NR1; bA350014.9 ;

Immunogen:

Human NDOR1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A