

NDUFA3 Antibody



PACO10776

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NDUFA3

Uniprot

O95167

Synonyms:

NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 3, 9kDa; NDUFA3; B9 ;

Immunogen:

Human NDUFA3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A