

# NDUFA5 Antibody



PACO10778

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NDUFA5

**Uniprot**

Q16718

**Synonyms:**

NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 5, 13kDa; NDUFA5; B13; CI-13KD-B; DKFZp781K1356; FLJ12147; NUFM; UQOR13 ;

**Immunogen:**

Human NDUFA5.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A