

NDUFA8 Antibody



PACO10780

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NDUFA8

Uniprot

P51970

Synonyms:

NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 8, 19kDa; NDUFA8; CI-19KD; CI-PGIV; MGC793; PGIV ;

Immunogen:

Human NDUFA8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A