

NDUFA9 Antibody



PACO10781

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NDUFA9

Uniprot

Q16795

Synonyms:

NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 9, 39kDa; NDUFA9; MGC111043; NDUFS2L; SDR22E1 ;

Immunogen:

Human NDUFA9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A