

NDUFB1 Antibody



PACO10784

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NDUFB1

Uniprot

O75438

Synonyms:

NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 1, 7kDa; NDUFB1; CI-SGDH; MNLL ;

Immunogen:

Human NDUFB1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A