

NDUFB2 Antibody



PACO10787

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NDUFB2

Uniprot

O95178

Synonyms:

NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 2, 8kDa; NDUFB2; AGGG; CI-AGGG; MGC70788 ;

Immunogen:

Human NDUFB2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A