

# NDUFB6 Antibody



PACO10789

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NDUFB6

**Uniprot**

O95139

**Synonyms:**

NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 6, 17kDa; NDUFB6; B17; CI; MGC13675 ;

**Immunogen:**

Human NDUFB6.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A