

NDUFB9 Antibody



PACO10792

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NDUFB9

Uniprot

Q9Y6M9

Synonyms:

NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 9, 22kDa; NDUFB9; B22; DKFZp566O173; FLJ22885; LYRM3; UQOR22 ;

Immunogen:

Human NDUFB9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A