

NECAB2 Antibody



PACO10802

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NECAB2

Uniprot

Q7Z6G3

Synonyms:

N-terminal EF-hand calcium binding protein 2; NECAB2; EFCBP2 ;

Immunogen:

Human NECAB2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A