

NECAP1 Antibody



PACO10804

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NECAP1

Uniprot

Q8NC96

Synonyms:

NECAP endocytosis associated 1; NECAP1; DKFZP566B183; MGC131900 ;

Immunogen:

Human NECAP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A