

NEDD8 Antibody



PACO10809

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NEDD8

Uniprot

Q15843

Synonyms:

neural precursor cell expressed, developmentally down-regulated 8; NEDD8; FLJ43224; MGC104393; MGC125896; MGC125897; Nedd-8 ;

Immunogen:

Human NEDD8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A