

# NEIL1 Antibody



PACO10812

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NEIL1

**Uniprot**

Q96FI4

**Synonyms:**

nei endonuclease VIII-like 1 (E. coli); NEIL1; FLJ22402; FPG1; NEI1; hFPG1 ;

**Immunogen:**

Human NEIL1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A