

NEK11 Antibody



PACO10814

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NEK11

Uniprot

Q8NG66

Synonyms:

NIMA (never in mitosis gene a)- related kinase 11; NEK11; FLJ23495 ;

Immunogen:

Human NEK11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A