

NEK2 Antibody



PACO10815

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NEK2

Uniprot

P51955

Synonyms:

NIMA (never in mitosis gene a)-related kinase 2; NEK2; HsPK21; NEK2A; NLK1 ;

Immunogen:

Human NEK2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A