

NEK9 Antibody



PACO10817

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NEK9

Uniprot

Q8TD19

Synonyms:

NIMA (never in mitosis gene a)- related kinase 9; NEK9; DKFZp434D0935; MGC138306; MGC16714; NERCC; NERCC1; Nek8 ;

Immunogen:

Human NEK9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A