

NELL1 Antibody



PACO10819

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NELL1

Uniprot

Q92832

Synonyms:

NEL-like 1 (chicken); NELL1; FLJ45906; IDH3GL; NRP1 ;

Immunogen:

Human NELL1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A