

NET1 Antibody



PACO10821

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NET1

Uniprot

Q7Z628

Synonyms:

596aa,68kd;

Immunogen:

Human NET1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

Product Images

N/A

N/A