

NETO2 Antibody



PACO10822

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NETO2

Uniprot

Q8NC67

Synonyms:

neuropilin (NRP) and tolloid (TLL)-like 2; NETO2; FLJ10430; FLJ14724; FLJ90456; NEOT2 ;

Immunogen:

Human NETO2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A