

NEURL1 Antibody



PACO10823

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NEURL1

Uniprot

O76050

Synonyms:

neuralized homolog (Drosophila); NEURL; NEURL1; RNF67; bA416N2.1; neu; neu-1 ;

Immunogen:

Human NEURL

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A