

NEURL2 Antibody



PACO10824

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NEURL2

Uniprot

Q9BR09

Synonyms:

neuralized homolog 2 (Drosophila); NEURL2; C20orf163; MGC125934; MGC125935; OZZ; Ozz-E3 ;

Immunogen:

Human NEURL2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A