

# NF2 Antibody



PACO10826

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NF2

**Uniprot**

P35240

**Synonyms:**

neurofibromin 2 (merlin); NF2; ACN; BANF; SCH ; Merlin; Schwannomin;  
Schwannomerlin; Neurofibromin-2;

**Immunogen:**

Human NF2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A