

NFE2L1 Antibody



PACO10829

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NFE2L1

Uniprot

Q14494

Synonyms:

nuclear factor (erythroid-derived 2)-like 1; NFE2L1; FLJ00380; LCR-F1; NRF1; TCF11 ;

Immunogen:

Human NFE2L1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A