

# NFKB1 Antibody



PACO10834

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NFKB1

**Uniprot**

P19838

**Synonyms:**

nuclear factor of  $\kappa$ ; light polypeptide gene enhancer in B-cells 1; NFKB1; DKFZp686C01211; EBP-1; KBF1; MGC54151; NF-kappa-B; NFKB-p105; NFKB-p50; p105; p50 ;

**Immunogen:**

Human NFKB1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20 $\text{deg}$ ; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A