

PACO10841

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NFYC

Uniprot

Q13952

Synonyms:

nuclear transcription factor Y, gamma; NFYC; CBF-C; CBFC; DKFZp667G242; FLJ45775; H1TF2A; HAP5; HSM; NF-YC ;

Immunogen:

Human NFYC.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A