

NKRF Antibody



PACO10863

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NKRF

Uniprot

O15226

Synonyms:

NFKB repressing factor; NKRF; ITBA4; NRF ;

Immunogen:

Human NKRF.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A