

NLN Antibody



PACO10868

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NLN

Uniprot

Q9BYT8

Synonyms:

neurolysin (metallopeptidase M3 family); NLN; AGTBP; DKFZp564F123; KIAA1226 ;

Immunogen:

Human NLN.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A