

# NLRP1 Antibody



PACO10869

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NLRP1

**Uniprot**

Q9C000

**Synonyms:**

NLR family, pyrin domain containing 1; NLRP1; CARD7; CLR17.1; DEFCAP; DEFCAP-L/S; DKFZp586O1822; KIAA0926; NAC; NALP1; PP1044; SLEV1; VAMAS1 ;

**Immunogen:**

Human NLRP1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A