

NME5 Antibody



PACO10876

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NME5

Uniprot

P56597

Synonyms:

non-metastatic cells 5, protein expressed in (nucleoside-diphosphate kinase); NME5; NM23-H5; NM23H5 ;

Immunogen:

Human NME5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A