

NMUR1 Antibody



PACO10880

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NMUR1

Uniprot

Q9HB89

Synonyms:

neuromedin U receptor 1; NMUR1; (FM-3); FM-3; FM3; GPC-R; GPR66; NMU1R ;

Immunogen:

Human NMUR1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A