

NOB1 Antibody



PACO10883

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

NOB1

Uniprot

Q9ULX3

Synonyms:

NIN1/RPN12 binding protein 1 homolog (S. cerevisiae); NOB1; ART-4; MST158; MSTP158; NOB1P; PSMD8BP1 ;

Immunogen:

Human NOB1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A