

# NOC4L Antibody



PACO10884

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NOC4L

**Uniprot**

Q9BVI4

**Synonyms:**

nucleolar complex associated 4 homolog (S. cerevisiae); NOC4L; MGC3162; NOC4 ;

**Immunogen:**

Human NOC4L.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A