

NOL6 Antibody



PACO10889

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NOL6

Uniprot

Q9H6R4

Synonyms:

nucleolar protein family 6 (RNA-associated); NOL6; FLJ21959; MGC14896; MGC14921; MGC20838; NRAP; UTP22; bA311H10.1 ;

Immunogen:

Human NOL6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A