

NOL8 Antibody



PACO10890

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NOL8

Uniprot

Q76FK4

Synonyms:

nucleolar protein 8; NOL8; C9orf34; DKFZp686P12242; FLJ20736; Nop132; bA62C3.3; bA62C3.4 ;

Immunogen:

Human NOL8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A