

NOMO1 Antibody



PACO10891

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NOMO1

Uniprot

Q15155

Synonyms:

NODAL modulator 1; NOMO1; Nomo; PM5 ;

Immunogen:

Human NOMO1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A