

NOS1 Antibody



PACO10897

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NOS1

Uniprot

P29475

Synonyms:

nitric oxide synthase 1 (neuronal); NOS1; IHPS1; NOS; nNOS ; NOS1

Immunogen:

Human nNOS.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A