

NOSTRIN Antibody



PACO10898

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NOSTRIN

Uniprot

Q8IV19

Synonyms:

nitric oxide synthase trafficker; NOSTRIN; DaIP2; MGC20702 ;

Immunogen:

Human NOSTRIN.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A