

PNP Antibody



PACO10903

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PNP

Uniprot

P00491

Synonyms:

nucleoside phosphorylase; NP; FLJ94043; FLJ97288; FLJ97312; MGC117396; MGC125915; MGC125916; PNP; PRO1837; PUNP ;

Immunogen:

Human NP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A