

NPAS1 Antibody



PACO10904

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NPAS1

Uniprot

Q99742

Synonyms:

neuronal PAS domain protein 1; NPAS1; MOP5; PASD5; bHLHe11 ;

Immunogen:

Human NPAS1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A