

NPC2 Antibody



PACO10906

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NPC2

Uniprot

P61916

Synonyms:

Niemann-Pick disease, type C2; NPC2; HE1; MGC1333; NP-C2 ;

Immunogen:

Human NPC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A