

# NQO2 Antibody



PACO10916

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NQO2

**Uniprot**

P16083

**Synonyms:**

NAD(P)H dehydrogenase, quinone 2; NQO2; DHQV; DIA6; NMOR2; QR2 ;

**Immunogen:**

Human NQO2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A