

NR1I2 Antibody



PACO10919

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NR1I2

Uniprot

O75469

Synonyms:

nuclear receptor subfamily 1, group I, member 2; NR1I2; BXR; ONR1; PAR; PAR1; PAR2; PARq; PRR; PXR; SAR; SXR ;

Immunogen:

Human NR1I2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A