

# NR113 Antibody



PACO10920

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NR113

**Uniprot**

Q14994

**Synonyms:**

nuclear receptor subfamily 1, group I, member 3; NR113; CAR; CAR1; MB67; MGC150433; MGC97144; MGC97209 ; CAR-BETA; Constitutive androstane receptor;

**Immunogen:**

Human NR113.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A