

# NR2C1 Antibody



PACO10921

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NR2C1

**Uniprot**

P13056

**Synonyms:**

nuclear receptor subfamily 2, group C, member 1; NR2C1; TR2 ;

**Immunogen:**

Human NR2C1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A