

NR2E3 Antibody



PACO10922

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NR2E3

Uniprot

Q9Y5X4

Synonyms:

nuclear receptor subfamily 2, group E, member 3; NR2E3; ESCS; MGC49976; PNR; RNR; RP37; rd7 ;

Immunogen:

Human NR2E3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A